



NORTH AMERICAN SCHOOL OF  
**PODOLOGY**



# Implements and Bits

For Foot Care Professionals

[www.NorthAmericanSchoolofPodology.com](http://www.NorthAmericanSchoolofPodology.com)

1-855-622-NASP (6277)

# Implements

Custom manufactured to meet the high quality standards of the North American School of Podology, the NASP Pedicure Implements are crafted with German precision from world renowned Solingen Stainless Steel. NASP is dedicated to providing individually handcrafted, highest quality implements known for precision, performance, longevity, reliability and uncompromising quality. Implements are meticulously handcrafted in a modern state of the art manufacturing facility using the best technologies and development capabilities. Implement experts inspect every implement using a 31 point check list guaranteeing quality and durability.

NASP implements are podiatry quality (surgical grade) without the podiatry price. With over 10 pieces to choose from, the NASP implements range from nail and cuticle nippers, probes, excavators, callus shavers and files. With a Lifetime Warranty\* against manufacturers' defects NASP Pedicure implements are a solid investment for nail technicians and salons.

NASP implements are made of podiatry grade cobalt stainless steel and can be chemically sterilized, autoclaved, or placed in dry a heat oven without damage. However, proper care and storage is imperative. Implements should be properly cleaned, disinfected or sterilized, dried and stored in an appropriate manner. NASP does not recommend submersing tools in disinfectants and chemosterilants that do not contain rust inhibitors (AHP not recommended).

## Cleaning/Disinfecting/Sterilizing

1. Clean implements with a soft brush and surfactant in cool water.
2. Rinse well and dry.
3. Submerge in disinfectant or chemosterilant for appropriate contact time, or place in autoclave or heat oven for appropriate contact time.
4. Rinse disinfectant or chemical sterilizer from tools in cool water.
5. Dry well especially in hinged areas.
6. Store in clean dry container.

***\*Every 2-3 months: oil hinges with oil for implements or surgical milk.***



# Nail and Skin Nippers

---



**NSP-209**  
**Toenail Clipper**  
Concave Jaw 14 cm

Controlled cutting of Toe Nails, excellent for cutting thick nails



**NSP-214**  
**Medium Flat Edge Nail Nipper**  
Flat edge jaw 15cm

Controlled cutting of sides , edges and spiculas



**NSP-117**  
**Small Flat Edge Nail Scissor**  
Flat edge jaw 10cm

Controlled trimming of cuticles, and removal of spiculas  
Small and easy to use for tight sidewalls

# Foot and Nail Files

---



**FLXSSDS-C/F**  
**Footlogix Double Sided Stainless Steel File**

Professional long lasting quality stainless steel, diamond system for multi directional filing without shredding the skin

# Excavators



**NSP-358  
Single Spoon**

Used for cleaning debris from nail groove and under free edge  
Small, smooth head for ease of use in tight areas and easy clean up



**NSP-410  
Double Spoon**



**NSP-311  
Corn Probe**

Precision tip for non-invasive corn reduction



**NSP-364  
Ingrown Nail File**

Narrow head for invasive removal of sharp edges and spiculae in the nail groove



**NSP-236  
Callus Plane**

Intricate head position for easy reduction of callus build-up

# NASP Implement Kit



**NSP-NASP-KIT-USA**



**NSP-NASP-KIT**

**Implements can be purchased separately or as a kit**

# Bits

---

NASP Bits are manufactured in Germany and are intended for medical pedicures, podiatry, and podology. All bits are constructed based on European Medical Device standards (Germany) using high quality stainless steel for the shanks, quality diamond, carbide and corundum for heads. Bits are meticulously handcrafted using the latest technologies and state of the art development capabilities for precision, performance, longevity, and uncompromising quality. The use of stainless steel ensures an optimal anti-rust protection for shanks. Diamond grit is galvanic bonded for long lasting quality.

All bits are surgical grade and safe for ultrasonic baths, heat ovens, autoclaves, and chemical disinfectants

## **Safety Recommendations:**

- Store in dry clean environment
- Insert properly into machine
- Make sure bits reach required speed before applying to skin or nail
- Only use maximum rotational speed recommended for each bit
- Bits should only be used in conformity with their intended use
- Wear personal protection equipment such as masks, safety goggles where required
- Do not use damaged bits – discard and replace immediately

## **Hygiene recommendations:**

- Clean debris with soft brush and surfactant under running water
- Can be used in ultrasonic bath
- Should not come in contact with one another during disinfecting or sterilizing
- AHP not recommended
- Make sure bits are dried with soft towel
- Sterilization in autoclave or heat oven recommended
- Proper storage imperative – use stand to keep implement in upright position for optimum bit care

NASP offers several options in bits from callus and nails reducers, to nail polishers and corn reducers. All bits are medical grade with affordable pricing for the pedicurist

# Callus Reducing Bits



**NSP-DT4880 120 HP**  
**Diatwister Mega-Coarse**



**NSP-DT5880 155 HP**  
**Diatwister Super-Coarse**

Initial cooling for working on callus  
Stainless steel shank with galvanized diamond heads  
Precise concentricity  
Use wet or dry



**NSP-5893 050**  
**Diamond Super**



**NSP-6894 063**  
**Diamond Coarse**



**NSP-8894 060**  
**Diamond Fine**

Used for precise callus reduction and smoothing along cracks and fissures  
Stainless steel shank with galvanized diamond heads  
Precise concentricity  
Use wet or dry



**NSP-HT6880 150**  
**Coarse Hybrid Grit**



**NSP-HT880 115**  
**Medium Hybrid Grit**



**NSP-HT8880 115**  
**Fine Hybrid Grit**

The Hybrid TWISTER's rate of removal will astound you. Two exceedingly effective materials, diamond and ceramic, combine to yield a progressive new hybrid grit. Diamond abrasive for a longer service life and higher cost-efficiency; Ceramic abrasive for a higher rate of removal and smooth surface

# Nail Reducers



**NSP-T429**  
**Speed Cutter**

With revolutionary cuts and safe-end working part, the new golden T-SPEED series of carbide cutters with Triple titanium coating produce optimum results quickly without heat.

The carbide cutter range SNC was specially developed for the treatment of thickened natural nails and artificial nail material for prosthetic corrections. A greater abrasion performance, low vibration, smooth surfaces, long service life and a safer handling



**NSP-854 050 HP**  
**Diamond Nail Reducer Large**

**NSP-854 KR 040**  
**Diamond Nail Reducer Medium**

**NSP-854 KR 033**  
**Diamond Nail Reducer Small**

Reduces mildly thick nails, precision callus work around nails and groove area.

Stainless steel shank, galvanized diamond head



**NSP-85RS 050 HP**  
**Steel Cut Bud**

Reduces thick gryphotic nails and hyperkeratoses

Coarse straight cut with cross cut  
Stainless steel shank and head

# Nail Sharpers and Smoothers



**NSP-840T 050**  
**Nail Shaper Coarse Top**  
**Grip**

**NSP-8840T 050**  
**Nail Shaper Fine Top**  
**Grip**

Shortening nails – safety  
bevelled at edge minimizing  
injury

Stainless Steel and Galvanized  
Diamond



**NSP-8840S 055**  
**Nail Plate Fine Side Grip**

Smoothly grinds top of plate,  
gentle filing right to nail wall

Bevelled edge minimizing risk  
of injury

Stainless Steel and Galvanized  
Diamond



**NSP-524 055 HP**  
**Abrasive White Small**  
**Nail Plate**

High grade corundum head, stainless steel shank  
Used for nail braces,  
Smoothing nail plate and shaping free edge



**NSP-524 100**  
**Abrasive White**  
**Large Nail Plate**

# Nail Plate Cleaners & Polishers



**NSP-648020 HP**  
**Abrasive Pink**

High grade corundum  
for work on nails and  
cuticle



**NSP-9628 060 HP**  
**Poliline Pol Anthraz**

Silicone bonded  
Polishes nail plate



**NSP-9624 060 HP**  
**Poliline Pol Khaki**

Silicone bonded  
Shines nail plate

# Specialty Bits



**NSP-8850 023 HP  
Diamond Fine for  
Nail Groove**

Fine diamond bit for reduction of callus in nail groove. Stainless steel shank, galvanized diamond head



**NSP-8863 012  
Diamond Fine  
Point for Nail  
Groove**



**NSP-39RS 009 HP  
Ingrown Nail File**

Intricate head for reduction of spicula in nail groove. Stainless steel shank and head



**NSP-11RS 018 HP  
Corn Reducer Small**

**NSP-11RS 009 HP  
Corn Reducer Large**

Precision round head for safe reduction of corns

Comes in 2 sizes for intricate reduction of corns



**NSP-15XSM 016 HP  
Round Gold Corn  
Reducer**

Precision round head for safe reduction of corns. Triple titanium for less heat

Gentle removal of corns and perigungal callus  
Slender neck allows better overview of area

# Bit Holder



**NSP-5764  
Display with Base**

Acrylic holder for storing bits.  
Fits over 15 bits  
Protects bits while holding them in an upright position



**NSP-5772  
Steri Safe Pro Bit Holder**

- Protects bits from damage while providing hygienic storage
- Holds up to 24 bits
- Made of quality durable plastic
- Lid locks to keep bits in place
- Ideal for mobile pedicurists



# B/S Brace Kit

Revolutionary bracing technique for the correction of involuted and ingrown toe nails. Made with fiberglass reinforced duroplast, the thin strips placed across the nail result in a gentle lifting, correcting the curvature of the nail. Pressure and pain is eliminated immediately.



## BS-010211 B/S Brace Classic Starter Set with Magnetic Application

### Kit Contains:

- 40 B/S-braces magnetic 2mm
- 10 each in the sizes 16, 18, 20, 22
- 5 B/S practice braces size 18
- 1 B/S Basic adhesive 5 g
- 1 B/S cleaner 25 ml
- 1 B/S Activator 10 ml
- 3 B/S magnetic applicators (all 3 sizes)
- 1 User manual
- 1 Counter Display
- 1 Informational DVD



## BS-104011 B/S Brace Classic+ Starter Set

### Kit Contains:

- 40 B/S-braces magnetic 4mm
- 10 each in the sizes 16, 18, 20, 22
- 5 B/S practice braces size 18
- 1 B/S Basic adhesive 5 g
- 1 B/S cleaner 25 ml
- 1 B/S Activator 10 ml
- 3 B/S magnetic applicators (all 3 sizes)
- 1 User manual
- 1 Counter Display
- 1 Informational DVD

All components can be purchased in kits or separately



### B/S Brace Clean-Ex 25 ml

BS Cleaner Clean EX is used to prepare the nail, improve adhesion and remove contaminants from the nail plate.



### B/S Brace Brush/Glue 5g

BS Basic-Glue is especially designed for use with the BS Brace. Simple application, quick adhesion and lasting wear.



### B/S Brace Magnetic Applicators

A magnet holds the BS Brace in place on the applicator, which makes it easier to apply the brace to the nail.



### B/S Brace Classic Aktivator

The BS Activator is used to activate the BS Glue and helps BS Brace last longer. Easy to use and quick to dry.

# Rechargeable E-File



## NSP-5764

- 30,000 RPM
- Rechargeable
- Long lasting charge (up to 8 hours)
- Variable speed control
- Light weight hand piece
- Quiet operation, no vibration
- Forward and reverse rotation

# Vacuum Suction E-File



## NSP-EFILEVAC

- Computerized control panel
- Programmable memory options for rotation and speed variables
- Slimmer hand piece with on/off touch control
- Reduced noise
- Reduced vibration
- Increased RPM 30,000
- Increased vacuum suction
- New improved cloth dust bag
- Full Dust Bag detection indicator



## NSP-0810 Micro Dust Bag for E-File

# Other Products

---



## **NSP-KANSFTWANDS**

### **Kansa Foot Wand**

Ancient Ayurveda massage device – wood handle, ball made from ZINC, COPPER & triple titanium ball “pulls” negative energy away from the sole of the foot, acidity is removed and foot chakras are balanced



## **NSP-PEDI BOWL LARGE**

### **Stainless Steel Foot Bowl Large**

Made from quality stainless steel, large for comfortably holding larger feet (16”x7”). Liners can be purchased to fit tubs



## **NSP-NASP # 953**

### **Instrument Sterilizing Tray with Cover**

Implement and bit disinfecting tray made from stainless steel. Removable tray with handles, lid fits snugly. Sized to hold foot files, one set of implements and bits